

STAIRWELL BELAY SYSTEM

The Stairwell Belay System offers complete control and safety for vertical evacuation. Controls the descent of the Evacuation Slyde/Bariatric Evacuation Slyde.



The Stairwell Belay System

For use with the Bariatric Evacuation Slyde and the Evacuation Slyde.



belluscura™

- **2 Stainless Steel Highly Durable Anchor Carabiners**

Connects the stairwell belay device and stairwell belay rope to a firm anchor point and to the Evacuation Slyde.

- **Stairwell Belay Rope**

30' of rope allows controlled descent down two flights of stairs.

- **Anchor Point Runner/Sling**

Wraps around a firm anchor point and connects a second time to the anchor carabiner.

- **Stairwell Belay Device**

Provides friction to control the descent. Minimizes the risk of injury during descent.

Works well with the
**STAIRWELL
BELAY BOX**

see other side for details

STAIRWELL BELAY BOX



Provides a stable wall anchor and storage for Stairwell Belay System.



Stairwell Belay Box
For use with the Bariatric Evacuation Slyde and the Stairwell Belay System.

(Evacuation Slyde and Stairwell Belay System pictured above)



belluscura™

- **Instant Storage**

Convenient instant storage of the Stairwell Belay System.

- **Locking Mechanism**

Keeps the Stairwell Belay System secure at all times. Includes key.

- **Stainless Steel Anchor Instant Access Point**

Provides a permanent, reliable, tamper-resistant anchor bolt access point in seconds for expedited evacuation.

- **Peace of Mind**

Reduces the risk of injury when descending the stairs.

ORDER INFORMATION

email: customerservice@belluscura.com

tel: +1 214-396-8313

toll free: +1 844-883-0940

Belluscura LLC

5504 Democracy Drive, Suite 200

Plano, TX 75024

United States

tel: +1 214-396-8313

Belluscura plc

82 St John Street, Farringdon

London EC1M 4JN

United Kingdom

tel: +44 (0) 1865 261506

toll free: +1 844-883-0940